



Your partner in safety and quality



BIOHOLE NEEDLE SYSTEM

Special designed needle for the
buttonhole puncture method



BioHole Needle System

For a painless arterial-venous insertion

The BioHole Needle System forms an excellent combination with the buttonhole puncture method. After a fix puncture tunnel is created, you can use the BioHole Needle System, which is less sharp than a standard AVF needle.

Buttonhole puncture may cause fewer complications, compared to using different needle insertion sites for each haemodialysis treatment. Therefore, the BioHole Needle System is recommended for home dialysis and self cannulation Renal Centers.

Using the BioHole Needle System in combination with the buttonhole puncture method results in :

- A reduction of bleeding time
- Minimum penetration force provides less pain to patients
- An extension of the life of the vascular access
- A reduction of mispunctures

Product Code	Needle		Back-Eye	Tube Length (mm)	Clamp Colour	Wing Colour
	Size (Gauge)	Length (mm-Inch)				
AV D 1520 HTC 15R	15	20 – ¾"	X	150	RED	
AV D 1520 HTC 30R	15	20 – ¾"	X	300	RED	
AV D 1520 TC 15B	15	20 – ¾"		150	BLUE	
AV D 1520 TC 30B	15	20 – ¾"		300	BLUE	
AV D 1525 HTC 15R	15	25 – 1"	X	150	RED	
AV D 1525 HTC 30R	15	25 – 1"	X	300	RED	
AV D 1525 TC 15B*	15	25 – 1"		150	BLUE	
AV D 1525 TC 30B*	15	25 – 1"		300	BLUE	
AV D 1620 HTC 15R*	16	20 – ¾"	X	150	RED	
AV D 1620 HTC 30R*	16	20 – ¾"	X	300	RED	
AV D 1620 TC 15B*	16	20 – ¾"		150	BLUE	
AV D 1620 TC 30B*	16	20 – ¾"		300	BLUE	
AV D 1625 HTC 15R*	16	25 – 1"	X	150	RED	
AV D 1625 HTC 30R*	16	25 – 1"	X	300	RED	
AV D 1625 TC 15B*	16	25 – 1"		150	BLUE	
AV D 1625 TC 30B*	16	25 – 1"		300	BLUE	
AV D 1725 HTC 15R*	17	25 – 1"	X	150	RED	
AV D 1725 HTC 30R*	17	25 – 1"	X	300	RED	
AV D 1725 TC 15B*	17	25 – 1"		150	BLUE	
AV D 1725 TC 30B*	17	25 – 1"		300	BLUE	

All needles have a coloured turnable wing and clamp
 Minimum order quantity for export of standard items 500 pcs
 * Contact your local office for order quantities and availability

NIPRO EUROPE N.V.
 European Headquarters
 Weihoek 3 H
 B-1930 Zaventem – Belgium
 Phone: +32 (0) 2 725 55 33
 Fax: +32 (0) 2 725 70 41
 e-mail: info@nipro-europe.be
 website: www.nipro-europe.com

NIPRO EUROPE N.V.
 Benelux Branch
 Weihoek 3H
 B-1930 Zaventem - Belgium
 Phone: +32 (0) 2 725 55 33
 Fax: +32 (0) 2 725 70 41

NIPRO EUROPE N.V.
 Italy Branch
 Centro Direzionale Milanofiori
 Strada 1 - Palazzo F/1
 I-20090 Assago (Milan) – Italy
 Phone: +39 02 57 50 00 57
 Fax: +39 02 57 51 81 11

NIPRO EUROPE N.V.
 Moscow Branch
 12 Krasnopresnenskaya nab.
 Office nbr 731 - 7th floor - entrance nr 6
 World Trade Center
 RF-123610 Moscow - Russian Federation
 Phone: +7 (8) 095 258 1366
 Fax: +7 (8) 095 258 1367

NIPRO EUROPE N.V.
 Spain Branch
 Avda. Castilla, 2 Edificio Berlin
 SP-28830 San Fernando de Henares, Madrid – Spain
 Phone: +34 91 675 18 34
 Fax: +34 91 675 19 17

NIPRO EUROPE N.V.
 Turkey Liaison Office
 Iskele Cad. Nil Apt. No. 17 - 3rd Floor - Flat 301
 34728 Caddebostan – Istanbul - Turkey
 Phone: +90 (0) 216 478 58 60
 Fax: +90 (0) 216 478 58 62

NIPRO EUROPE N.V.
 U.K. Branch
 Ground Floor – Unit 7
 The Forum
 Minerva Business Park Lynchwood
 GB-Peterborough PE2 6FT – United Kingdom
 Phone : +44 (0) 173 323 4300
 Fax : +44 (0) 173 323 4310

